

Hoechst Staining protocol for cells

Hoechst 33528 stain: Sigma 14530

Make 1 mg/ml solution in ddH₂O

Store protected from light at 4 deg C

On day of assay, dilute to 20 ug/mL in buffer TNE (10 mM Tris, 1 mM EDTA, 2M NaCl, pH 7.4, stored at RT)

After assay, remove media from cell plate by overturning over paper towels

Freeze plate at -80deg C

On day of assay, remove plate and add 100 ul ddH₂O per well; incubate at RT for 1 hour

Freeze plate for >20 mins at -80 deg C

Remove plate and allow to thaw to room temperature

Add 100 ul of 20ug/mL Hoechst in TNE buffer for final concentration of 10 ug/mL

Measure fluorescence: Excite 350 nm, measure emission at 460 nm.

For questions, contact Harini Sampath, sampath.harini@gmail.com